



Anti-Hexestrol polyclonal antibody (DPAB-DC4348)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	ethylstilbestrol 43%, lbestrol-Glucuronide 85%, Dienestrol 20%
Immunogen	Hexestrol conjugated with HSA.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Hexestrol is a carcinogenic synthetic estrogen that inhibits microtubule polymerization and the formation of ribbon structures. Hexestrol is also an inhibitor of lipid peroxidation.
Keywords	Hexestrol