



Anti-Histone H1 polyclonal antibody (DPAB-DC4355)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Intact calf histone H1 complexed with RNA
Source/Host	Sheep
Species Reactivity	Bovine
Conjugate	Unconjugated
Applications	WB (Recombinant protein), IF, IP,
Format	Liquid
Size	250 µg
Buffer	In PBS, pH 7.4
Preservative	None
Storage	Store at -20°C for short term. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Five major families of histones exist: H1/H5, H2A, H2B, H3 and H4. Histones H2A, H2B, H3 and H4 are known as the core histones, while histones H1 and H5 are known as the linker histones.
Keywords	Histone H1