



Anti-Histone H3 polyclonal antibody (DPAB-DC4363)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) containing the dimethylated lysine 9 (or [K9me2]) of human histone H3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ChIP, WB (Cell lysate), ELISA, Dot,
Format	Liquid
Size	50 µg
Buffer	In PBS (0.05% sodium azide, 0.05% proclin 300)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Five major families of histones exist: H1/H5, H2A, H2B, H3 and H4. Histones H2A, H2B, H3 and H4 are known as the core histones, while histones H1 and H5 are known as the linker histones.
--------------	--

Keywords

H3K9me2; Histone H3 (K9me2)
