



Anti-HMGA2-P3 polyclonal antibody (DPAB-DC4376)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to HMGA2-P3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, IHC,
Format	Liquid
Size	100 µg
Buffer	In 15 mM KPO4 buffer, 0.85% NaCl, pH 7.4 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. Do not freeze.

BACKGROUND

Keywords	HMGA2-P3
----------	----------

GENE INFORMATION

Entrez Gene ID [84017](#)

UniProt ID [O88791](#)
