



Anti-HSP101 (N-terminal) polyclonal antibody (DPAB-DC4381)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Detects ~101KDa.
Immunogen	Recombinant protein corresponding to N-terminus of Arabidopsis thaliana HSP101.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Plants, Arabidopsis thaliana
Conjugate	Unconjugated
Applications	WB (Tissue lysate),
Format	Liquid
Size	100 µl
Buffer	In antiserum
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	HSP101
----------	--------