



Anti-HSP90 (full length) polyclonal antibody (DPAB-DC4385)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody recognizes eukaryotic Hsp90alpha/beta.
Immunogen	Full length Hsp90 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µl
Buffer	In serum
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Hsp90 (heat shock protein 90) is a chaperone protein that assists other proteins to fold properly, stabilizes proteins against heat stress, and aids in protein degradation. It also stabilizes a number of proteins required for tumor growth, which is why Hs
--------------	---

