



Anti-HHV type 2 ICP34.5 Polyclonal antibody (DPAB-DC4390)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	HHV type 2 ICP34.5
Immunogen	A synthetic peptide corresponding to N-terminus of human herpesvirus 2 ICP34.5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HHV
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Herpesviridae is a large family of DNA viruses that cause diseases in animals, including humans. The members of this family are also known as herpesviruses. The family name is derived from the Greek word herpein ("to creep"), referring to the latent, recu
--------------	---

Keywords

HHV-2 ICP34.5; Human herpesvirus type 2 ICP34.5
