



Anti-HIV Glycoprotein 120 Polyclonal antibody (DPAB-DC4396)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	HIV Glycoprotein 120
Immunogen	A synthetic peptide corresponding to the internal region 15 amino acids of HIV gp120.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HIV
Conjugate	Unconjugated
Applications	WB (Recombinant protein),
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The human immunodeficiency virus (HIV) is a lentivirus (a subgroup of retrovirus) that causes the acquired immunodeficiency syndrome (AIDS), a condition in humans in which progressive failure of the immune system allows life-threatening opportunistic infe
--------------	---

Keywords

HIV gp120; Human immunodeficiency virus gp120
