



Anti-HIV P24 Polyclonal antibody (DPAB-DC4397)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody recognizes both native and recombinant HIV-P24.
Target	HIV P24
Immunogen	Recombinant protein corresponding to human Immunodeficiency Virus (HIV) P24.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HIV
Conjugate	Unconjugated
Applications	ELISA,
Format	Lyophilized
Size	1 mg
Buffer	Lyophilized from 20 mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C. Aliquot after reconstitution to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The human immunodeficiency virus (HIV) is a lentivirus (a subgroup of retrovirus) that causes the acquired immunodeficiency syndrome (AIDS), a condition in humans in which progressive
--------------	---

failure of the immune system allows life-threatening opportunistic infe

Keywords

HIV P24; Human immunodeficiency virus P24
