



# Anti-HIV type 1 Gag p17 Polyclonal antibody (DPAB-DC4400)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	HIV type 1 Gag p17
Immunogen	Recombinant protein corresponding to human immunodeficiency virus type 1 Gag p17.
Source/Host	Pig
Species Reactivity	HIV
Conjugate	Unconjugated
Applications	WB (Cell lysate), ELISA,
Format	Liquid
Size	25 µl
Buffer	In buffer containing 0.09% sodium azide
Preservative	0.09% Sodium Azide
Storage	Store at 4°C.

## BACKGROUND

Introduction	HIV-1 is the most common and pathogenic strain of the virus. Scientists divide HIV-1 into a major group (Group M) and two or more minor groups. Each group is believed to represent an independent transmission of SIV into humans (but subtypes within a group
Keywords	HIV-1 Gag p17; Human immunodeficiency virus type 1 Gag p17