



Anti-HIV type 1 Polyclonal antibody (DPAB-DC4408)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody recognizes: P17, P24, P45, P53, P55, P64, P120 and gp160. Strong positive bands on P17, P24 and P64.
Target	HIV type 1
Immunogen	Native purified human immunodeficiency virus (HIV).
Source/Host	Goat
Species Reactivity	HIV
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Lyophilized
Size	1 mg
Buffer	Lyophilized from 20 mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C. Aliquot after reconstitution to avoid repeated freezing and thawing.

BACKGROUND

Introduction	HIV-1 is the most common and pathogenic strain of the virus. Scientists divide HIV-2 into a major group (Group M) and two or more minor groups. Each group is believed to represent an independent transmission of SIV into humans (but subtypes within a group
---------------------	---

Keywords

HIV-1; Human immunodeficiency virus type 1
