



Anti-HIV type 1 Reverse Transcriptase Polyclonal antibody (DPAB-DC4410)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	HIV type 1 Reverse Transcriptase
Immunogen	Recombinant protein corresponding to human immunodeficiency virus type 1 reverse transcriptase.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HIV
Conjugate	Unconjugated
Applications	WB (Cell lysate), ELISA,
Format	Liquid
Size	50 µl
Buffer	In buffer containing 0.09% sodium azide
Preservative	0.09% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction	HIV-1 is the most common and pathogenic strain of the virus. Scientists divide HIV-4 into a major group (Group M) and two or more minor groups. Each group is believed to represent an independent transmission of SIV into humans (but subtypes within a group
---------------------	---

Keywords

HIV-1 RT; Human immunodeficiency virus type 1 Reverse Transcriptase
