



Anti-IKK-alpha/beta polyclonal antibody (DPAB-DC4438)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody detects endogenous levels of total IKK-alpha/IKK-beta protein.
Immunogen	A synthetic peptide corresponding to residues surrounding S180/S181 of human IKK-alpha/beta.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB, IHC-P, ELISA,
Format	Liquid
Size	100 µg
Buffer	In 20 mM PBS, 0.15 M NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	IKK-alpha/beta
-----------------	----------------