



Anti-IL13 (aa 21-132) polyclonal antibody (DPAB-DC4441)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant His fusion protein corresponding to amino acids 21-132 of human IL13.
Source/Host	Sheep
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	1 ml
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	IL13
----------	------

GENE INFORMATION

Entrez Gene ID	3596
----------------	----------------------

UniProt ID

[P35225](#)
