



Anti-IL1F5 polyclonal antibody (DPAB-DC4442)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to 155 amino acids of bovine IL1F5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	IL1F5
-----------------	-------

GENE INFORMATION

Entrez Gene ID	54450
-----------------------	-----------------------

UniProt ID

[A2AIA9](#)
