



Anti-IL21 (aa 24-152) polyclonal antibody (DPAB-DC4443)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids 24-152 of bovine IL21.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine, Horse, Pig
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Format	Liquid
Size	100 µg
Buffer	In 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	IL21
----------	------

GENE INFORMATION

Entrez Gene ID [59067](#)

UniProt ID [Q9HBE4](#)
