



Anti-INA polyclonal antibody (DPAB-DC4445)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Purified bovine INA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine, Mammals
Conjugate	Unconjugated
Applications	WB, ICC, IF,
Format	Liquid
Size	100 µl
Buffer	In PBS (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	INA
----------	-----

GENE INFORMATION

Entrez Gene ID	9118
UniProt ID	Q16352