



# Anti-IAV H1N1 Polyclonal antibody (DPAB-DC4448)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Purified virions. Specific to H1N1 by IHA. Does not react with Influenza B, RSV, Parainfluenza 1-3 or Adenovirus. Does not react with HEp-2 cells. May react with chicken cellular proteins.
<b>Target</b>	IAV H1N1
<b>Immunogen</b>	Influenza A, strain A/USSR/90/77 (H1N1).
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	IAV
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P, IF,
<b>Format</b>	Liquid
<b>Size</b>	1 ml
<b>Buffer</b>	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Influenza A virus causes influenza in birds and some mammals, and is the only species of influenza virus A. Influenza virus A is a genus of the Orthomyxoviridae family of viruses. Strains
---------------------	---

of all subtypes of influenza A virus have been isolated from wild b

---

**Keywords**

H1N1 Virions; Influenza A H1N1 (Virions)

---