



Anti-IAV Polyclonal antibody (DPAB-DC4451)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Purified virions. Does not cross-react with influenza A virus, Parainfluenza virus (types 1-3), Respiratory Syncytial virus or Adenovirus. Negative against monkey kidney cells by indirect immunofluorescence.
Target	IAV
Immunogen	Influenza B, strain Yamagata.
Source/Host	Goat
Species Reactivity	IAV
Conjugate	Unconjugated
Applications	WB, IF, ELISA,
Format	Liquid
Size	100 µl
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Influenza A virus causes influenza in birds and some mammals, and is the only species of influenza virus A. Influenza virus A is a genus of the Orthomyxoviridae family of viruses. Strains of all subtypes of influenza A virus have been isolated from wild b
---------------------	---

Keywords

virions; Influenza B (virions)
