



# Anti-IAV Hemagglutinin Polyclonal antibody (DPAB-DC4452)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	IAV Hemagglutinin
<b>Immunogen</b>	A synthetic peptide (conjugated with KLH) corresponding to influenza hemagglutinin-HA. The sequence is YPYDVPDYA
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	IAV
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS, pH 7.2
<b>Preservative</b>	None
<b>Storage</b>	Store at 4°C.

## BACKGROUND

<b>Introduction</b>	Influenza A virus causes influenza in birds and some mammals, and is the only species of influenza virus A. Influenza virus A is a genus of the Orthomyxoviridae family of viruses. Strains of all subtypes of influenza A virus have been isolated from wild b
---------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**Keywords**

HA; Influenza hemagglutinin-HA

---