



Anti-S. frugiperda Cell Line SF21 Polyclonal antibody (DPAB-DC4455)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts with >40 SF21 proteins by Western blot.
Target	S. frugiperda Cell Line SF21
Immunogen	Native purified from the ovarian derived cell line of Spodoptera frugiperda.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. frugiperda
Conjugate	Unconjugated
Applications	WB, ELISA, Dot,
Format	Liquid
Size	1 ml
Buffer	In PBS
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Insect Cell Line SF21
----------	-----------------------