



# Anti-S. frugiperda Cell Line SF9 Polyclonal antibody (DPAB-DC4456)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Reacts with more than 40 SF9 proteins by Western blot.
Target	S. frugiperda Cell Line SF9
Immunogen	Native purified from the ovarian derived cell line of Spodoptera frugiperda.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. frugiperda
Conjugate	Unconjugated
Applications	WB, ELISA, Dot,
Format	Liquid
Size	1 ml
Buffer	In citrate-phosphate buffer, pH 7.0-7.4
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	Insect Cell Line SF9
----------	----------------------