



Anti-Keratin 12 polyclonal antibody (DPAB-DC4465)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to rabbit keratin 12.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rabbit
Conjugate	Unconjugated
Applications	IHC, IF,
Format	Liquid
Size	25 µg
Buffer	In PBS (0.1% proclin, 2.0% Block Ace)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Keratin 12 is a protein that in humans is encoded by the KRT12 gene. Keratin 12 is keratin found expressed in corneal epithelia. Mutations in the gene encoding this protein lead to Meesmann corneal dystrophy.
Keywords	Keratin 12

GENE INFORMATION

Entrez Gene ID [3859](#)

UniProt ID [Q99456](#)
