



Anti-Keratin (native) polyclonal antibody (DPAB-DC4466)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes keratin.
Immunogen	Native purified human Keratin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Birds, Human, Mammals
Conjugate	Unconjugated
Applications	IHC,
Format	Liquid
Size	250 µl
Buffer	In purified antiserum
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Keratin is a family of fibrous structural proteins. Keratin is the key structural material making up the outer layer of human skin. It is also the key structural component of hair and nails. Keratin monomers assemble into bundles to form intermediate fila
--------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Keywords

Keratin

GENE INFORMATION

Entrez Gene ID

[3852](#)

UniProt ID

[P13647](#)
