



Anti-Ketoprofen polyclonal antibody (DPAB-DC4467)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Ketoprofen conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Ketoprofen, (RS)2-(3-benzoylphenyl)-propionic acid (chemical formula C ₁₆ H ₁₄ O ₃) is one of the propionic acid class of nonsteroidal anti-inflammatory drugs (NSAID) with analgesic and antipyretic effects. It acts by inhibiting the body's production of prosta
Keywords	Ketoprofen