



# Anti-Ki67 antigen polyclonal antibody (DPAB-DC4468)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant protein corresponding to a 1086 bp Ki67 motif-containing cDNA fragment of Ki67 Antigen.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P,
<b>Format</b>	Lyophilized
<b>Size</b>	200 µl
<b>Buffer</b>	Lyophilized from PBS (1% BSA, 15 mM sodium azide)
<b>Preservative</b>	15mM Sodium Azide
<b>Storage</b>	Store at 4°C on dry atmosphere. After reconstitution with steril distilled water, Store at 4°C for at least two months. For long term storage store at -20°C or below. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

**Keywords** Ki67 Antigen

# GENE INFORMATION

Entrez Gene ID [4288](#)

---

UniProt ID [P46013](#)

---