



Anti-Klebsiella Polyclonal antibody (DPAB-DC4469)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts primarily with <i>K. pneumoniae</i> but cross-reacts with other <i>Klebsiella</i> species as well.
Target	<i>Klebsiella</i>
Immunogen	Native genus specific proteins from <i>Klebsiella pneumoniae</i> .
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Klebsiella</i>
Conjugate	Unconjugated
Applications	sELISA, Capture Antibody, Agg, Detection Antibody,
Format	Liquid
Size	1 ml
Buffer	In buffer containing 0.1% proclin 300
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction *Klebsiella* species are ubiquitous in nature. This is thought to be due to distinct sublineages

developing specific niche adaptations, with associated biochemical adaptations which make them better suited to a particular environment. They can be found in w

Keywords

Klebsiella; Klebsiella (Genus Specific)
