



Anti-KLH polyclonal antibody (DPAB-DC4471)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Purified KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA,
Format	Lyophilized
Size	50 µg
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at -20°C on dry atmosphere. After reconstitution with sterile water, store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Keyhole limpet hemocyanin (KLH) is a large, multisubunit, oxygen-carrying, metalloprotein that is found in the hemolymph of the giant keyhole limpet, <i>Megathura crenulata</i> , a species of keyhole limpet that lives off the coast of California, from Monterey Ba
Keywords	KLH