



Anti-KLH polyclonal antibody (DPAB-DC4473)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Keyhole Limpet Hemocyanin, KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, IHC, ELISA,
Format	Liquid
Size	1 mg
Buffer	In 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Keyhole limpet hemocyanin (KLH) is a large, multisubunit, oxygen-carrying, metalloprotein that is found in the hemolymph of the giant keyhole limpet, <i>Megathura crenulata</i> , a species of keyhole limpet that lives off the coast of California, from Monterey Ba
Keywords	KLH