



Anti-KLH polyclonal antibody (DPAB-DC4474)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A mixture of synthetic peptides (conjugated to KLH) corresponding to KLK5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB,
Format	Lyophilized
Size	200 µl
Buffer	Lyophilized from serum
Preservative	None
Storage	Store at -20°C on dry atmosphere. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Keyhole limpet hemocyanin (KLH) is a large, multisubunit, oxygen-carrying, metalloprotein that is found in the hemolymph of the giant keyhole limpet, Megathura crenulata, a species of keyhole limpet that lives off the coast of California, from Monterey Ba
Keywords	KLH; KLK5

GENE INFORMATION

Entrez Gene ID [25818](#)

UniProt ID [Q9Y337](#)
