



Anti-Laminin (native) polyclonal antibody (DPAB-DC4483)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	No immunoreactivity was obtained with human type I, III, IV and V collagen, whereas this antibody reacted with a distinct band of approximately 200 - 220 kD from a 8 M Urea extract from amion basement membrane.
Immunogen	Native purified mouse Laminin .
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mammals, Mouse, Rat
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr,
Format	Liquid
Size	50 µl
Buffer	In serum (sodium azide)
Preservative	See individual product datasheet
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Laminins are major proteins in the basal lamina (one of the layers of the basement membrane),
---------------------	---

a protein network foundation for most cells and organs. The laminins are an important and biologically active part of the basal lamina, influencing cell differen

Keywords	Laminin
-----------------	---------

GENE INFORMATION

Entrez Gene ID	3913
-----------------------	----------------------

UniProt ID	P55268
-------------------	------------------------
