



Anti-Lasalocid A polyclonal antibody (DPAB-DC4485)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Lasalocid A conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Lasalocid is able to make neutral complexes with monovalent and divalent cations and transport them through apolar phase (including lipid bilayer membranes). It can transport also big organic cations like dopamine.
Keywords	Lasalocid A