



Anti-L-aspartic acid polyclonal antibody (DPAB-DC4486)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	L-aspartic acid conjugated with carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	IHC, ICC, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Aspartic acid (abbreviated as D-AA, Asp, or D) is an α -amino acid with the chemical formula $\text{HOOCCH}(\text{NH}_2)\text{CH}_2\text{COOH}$. The carboxylate anion and salts of aspartic acid are known as aspartate. The L-isomer of aspartate is one of the 22 proteinogenic amino acids,
--------------	---

Keywords

L-aspartic acid
