



Anti-LDL (native protein) polyclonal antibody (DPAB-DC4492)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified human low density lipoprotein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, IP, EIA,
Format	Liquid
Size	100 µg, 150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Low-density lipoprotein (LDL) is one of the five major groups of lipoproteins, which, in order of molecular size, largest to smallest, are chylomicrons, very low-density lipoprotein (VLDL), intermediate-density lipoprotein (IDL), LDL, and high-density lip
Keywords	LDL