



Anti-L. biflexa Polyclonal antibody (DPAB-DC4493)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Leptospira biflexa, all antigens. Reactive in ELISA with L. biflexa, L. australis, L. canicola, L. grippotyphosa, L. icterohaemorrhagiae, L. pomona, L. serjoe. Antiserum has not been absorbed and may react with related microorganisms.
Target	L. biflexa
Immunogen	Native from Leptospira biflexa.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	L. biflexa
Conjugate	Unconjugated
Applications	IF, ELISA, EIA,
Format	Liquid
Size	1 ml
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords

L. biflexa; Legionella biflexa
