



Anti-L. pneumophila Polyclonal antibody (DPAB-DC4494)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts with L. pneumophila serogroups 1-12 in IFA. Recognizes all antigens of intact microorganism. Antiserum is not absorbed and may cross-react with related microorganisms.
Target	L. pneumophila
Immunogen	Native from Legionella pneumophila (ATCC #33152).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	L. pneumophila
Conjugate	Unconjugated
Applications	IF, EIA,
Format	Liquid
Size	100 µl
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction

Legionella pneumophila is a thin, aerobic, pleomorphic, flagellated, non-spore forming, Gram-negative bacterium of the genus Legionella. L. pneumophila is the primary human pathogenic bacterium in this group and is the causative agent of legionellosis or

Keywords

L. pneumophila; Legionella pneumophila
