



Anti-LOC440866 (aa 2-88) polyclonal antibody (DPAB-DC4503)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	LOC440866 (XP_496551, 2 a.a. ~ 88 a.a) partial recombinant protein with GST tag. The sequence is ATHGQTCARPMCIPPSYADLGKAARDIFNKGFGFLVKLDVKT KSCSGVEFSTSGSSNTD TGKVTGTLET KYKWCEYGLTFTEKWNT
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	LOC440866; LOC440866 (A01)
-----------------	----------------------------