



# Anti-MAM1 polyclonal antibody (DPAB-DC4509)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |  |
|--------------------|--|
| Immunogen          | A synthetic peptide corresponding to human MAM1.   |
| Isotype            | IgG  |
| Source/Host        | Rabbit   |
| Species Reactivity | Yeast  |
| Conjugate          | Unconjugated   |
| Applications       | WB,  |
| Format             | Liquid   |
| Size               | 100 µg   |
| Buffer             | In PBS (0.02% sodium azide)  |
| Preservative       | 0.02% Sodium Azide   |
| Storage            | Store at 4°C for short term. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing. |

## BACKGROUND

|          |      |
|----------|------|
| Keywords | MAM1 |
|----------|------|