



Anti-MAPK1/MAPK3 polyclonal antibody (DPAB-DC4510)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to 35 amino acids of MAPK1/MAPK3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine, Chicken, Frog, Fruit fly, Human, Mouse, Rat, Sheep
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µl
Buffer	In PBS
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords MAPK1/MAPK3

GENE INFORMATION

Entrez Gene ID [5594](#)

UniProt ID [P28482](#)
