



Anti-Melanophoren Stimulierendes Hormon polyclonal antibody (DPAB-DC4518)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody reacts with alpha-MSH (Melanophore Stimulating Hormone).
Immunogen	A synthetic peptide (conjugated with BSA) corresponding to Melanophoren stimulierendes hormon.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Pig, Rat
Conjugate	Unconjugated
Applications	IHC,
Format	Liquid
Size	250 µl
Buffer	In serum
Preservative	None
Storage	Store at 4°C.

BACKGROUND

Keywords	Melanophoren stimulierendes hormon; Melanophoren Stimulierendes Hormon
-----------------	--