



Anti-Melengestrol polyclonal antibody (DPAB-DC4521)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Melengestrol conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Melengestrol (INN, BAN) is a steroidal progestin and antineoplastic agent which was never marketed. An acylated derivative, melengestrol acetate, is used as a growth promoter in animals. While melengestrol is sometimes used as a synonym for melengestrol a
Keywords	Melengestrol