



Anti-Meloxicam polyclonal antibody (DPAB-DC4522)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Meloxicam conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Meloxicam is a nonsteroidal anti-inflammatory drug (NSAID) with analgesic and fever reducer effects. It is a derivative of oxicam, closely related to piroxicam, and falls in the enolic acid group of NSAIDs. It was developed by Boehringer-Ingelheim. Meloxi
Keywords	Meloxicam