



Anti-Mestanolone polyclonal antibody (DPAB-DC4525)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Mestanolone conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Mestanolone is the 17 α -methylated version of dihydrotestosterone (DHT). The systematic name of mestanolone is: 17 β -hydroxy-17 α -methylandroster-3-one. It is an orally bioavailable androgenic steroid that is highly androgenic while only slightly anabolic. It i
Keywords	Mestanolone