



Anti-S. mansoni Metabotropic glutamate receptor Polyclonal antibody (DPAB-DC4526)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	S. mansoni Metabotropic glutamate receptor
Immunogen	A synthetic peptide corresponding to C-terminus of Schistosoma mansoni metabotropic glutamate receptor.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. mansoni
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The metabotropic glutamate receptors, or mGluRs, are a type of glutamate receptor that are active through an indirect metabotropic process. They are members of the group C family of G-
---------------------	--

protein-coupled receptors, or GPCRs. Like all glutamate receptors, mGlu

Keywords

mGluR; metabotropic glutamate receptor
