



Anti-MGC42638 (aa 1-75) polyclonal antibody (DPAB-DC4539)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	MGC42638 (AAH40290, 1 a.a. ~ 75 a.a) partial recombinant protein with GST tag. The sequence is MALKLIHKEFLELARDPQPHCSAGPVWDDMLHWQATITRPNDSSYLGGVFFLKFPSDYLF KPPKIKFTNGIYHQR
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	MGC42638; MGC42638 (A01)
----------	--------------------------