



Anti-MLC1 (N-terminal) polyclonal antibody (DPAB-DC4540)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus of Alaska pollock MLC1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fish
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Alaska pollock or walleye pollock (<i>Theragra chalcogramma</i>) is a species of the cod family Gadidae. Alaska pollock is a semipelagic schooling fish widely distributed in the North Pacific with largest concentrations found in the eastern Bering Sea.
Keywords	MLC1

GENE INFORMATION

Entrez Gene ID [23209](#)

UniProt ID [Q15049](#)
