



Anti-MLL4 (internal region) polyclonal antibody (DPAB-DC4541)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A mixture of synthetic peptides (conjugated with KLH) that corresponding to internal region of mouse Mll4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB (Cell lysate), ELISA,
Format	Liquid
Size	100 µl
Buffer	In serum (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords Mll4

GENE INFORMATION

Entrez Gene ID [9757](#)

UniProt ID [Q9UMN6](#)
