



Anti-MMP-9 polyclonal antibody (DPAB-DC4542)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody does not react with other MMPs in IHC.
Immunogen	A synthetic peptide corresponding to MMP-9 (Gelatinase B).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, IHC,
Format	Liquid
Size	1 ml
Buffer	In PBS (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Gelatinase B (EC 3.4.24.35, 92-kDa gelatinase, matrix metalloproteinase 9, type V collagenase, 92-kDa type IV collagenase, macrophage gelatinase, 95 kDa type IV collagenase/gelatinase, collagenase IV, collagenase type IV, gelatinase MMP 9, MMP 9, type IV
---------------------	--

Keywords

MMP-9; MMP-9 (Gelatinase B)
