



Anti-Monensin polyclonal antibody (DPAB-DC4544)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Monensin conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Monensin is a polyether antibiotic isolated from <i>Streptomyces cinnamonensis</i> . It is widely used in ruminant animal feeds.
Keywords	Monensin