



Anti-M. morganii Polyclonal antibody (DPAB-DC4546)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	M. morganii
Immunogen	Morganella morganii total protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	M. morganii
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction Morganella morganii is a species of Gram-negative bacteria. M. morganii has a commensal

relationship within the intestinal tracts of humans, mammals, and reptiles as normal flora. Although *M. morganii* has a wide distribution, it is considered an uncommon

Keywords	<i>M. morganii</i> ; <i>Morganella morganii</i>
-----------------	---
